



КАРТА ФАЙЛОВ / TORQUE CARD





ProTaper Universal™
www.dentsplymaillefer.com

Файл / File		SP, rpm	TQ, Ncm
Sx	-	150-350 (250)	3,0
S1	1		
S2	1		1,0
F1 (20/.07)	1		1,5
F2 (25/.08)	1		2,0
F3 (30/.09)	1		
F4 (40/.06)	2		
F5 (50/.05)	2		
D1 (30/.09)	1		500* / (250 -300)**
D2 (25/.08)	2		
D3 (20/.07)	3	1,5	

*for obturations with gutta percha or Thermafil/
ProTaper® obturators

** for zinc oxide eugenol based obturations



ProFile®
www.dentsplymaillefer.com

Файл / File		SP, rpm	TQ, Ncm	
O.S.#6 (80/.08)	3	150-350 (300)	3,0	
O.S.#5 (60/.08)	3			
O.S.#4 (50/.07)	3			2,0
O.S.#3 (40/.06)	3			
O.S.#2 (30/.06)	3			
O.S.#1 (20/.05)	3			
40/.06	2		1,5	
35/.06	2			
30/.06	2			
25/.06	2			
20/.06	2			
15/.06	2			
90/.04	1		3,0	
60/.04	1			
45/.04	1		2,0	
40/.04	1		1,0	

Файл / File		SP, rpm	TQ, Ncm
35/.04	1	150-350 (300)	1,0
30/.04	1		0,6
25/.04	1		
20/.04	1		
15/.04	1		
40/.02	1		1,0
35/.02	1		0,6
30/.02	1		
25/.02	1		
20/.02	1		
15/.02	1		



GT® Series X™
www.dentsplymaillefer.com

Файл / File			SP, rpm	TQ, Ncm
20/.06		III	300	1,0
20/.04		II		
30/.08		III		2,0
30/.06		III		1,5
30/.04		II		1,0
40/.08		III		2,0
40/.06		III		1,5
40/.04		II		



PathFile™
www.dentsplymaillefer.com

Файл / File		SP, rpm	TQ, Ncm
013	1	300	max
016	1		
019	1		

Файл / File			SP, rpm	TQ, Ncm
INTRO FILE		-	250 - 350 (280)	1,4
.06/40		III		1,9
.06/35		III		1,6
.06/30		III		1,3
.06/25		III		0,9
.06/20		III		0,7
.06/15		III		0,4
.04/40		II		1,3
.04/35		II		1,0
.04/30		II		0,9
.04/25		II		0,6
.04/20		II		0,4
.04/15		II		0,3
.02/70		I		3,5
.02/60		I		2,5
.02/50		I		1,9

Файл / File			SP, rpm	TQ, Ncm
.02/45		I	250 - 350 (280)	1,6
.02/40		I		1,1
.02/35		I		0,7
.02/30		I		0,5
.02/25		I		0,4
.02/20		I		0,2
.02/15		I		

Файл / File			SP, rpm	TQ, Ncm
10/.04		I	250 - 350 (280)	1,2
15/.05		II		1,3
20/.06		III		2,1
25/.06		III		2,3
30/.05		II		1,2
35/.04		I		
40/.04		I		1,6
45/.04		I		
50/.04		I		2,0
60/.04		I		3,0
25/.07		III		2,0
30/.06		III		1,2
35/.06		III		1,0
40/.06		III		1,7
R15/.05				0,3
R25/.05				1,2



EasyShape™ System
www.kometdental.de

Файл / File		SP, rpm	TQ, Ncm
10.035		500 600	2,0 1,6
04.015		250-500	0,2
05.020			0,3
06.025			0,7
05.030			
04.035			
04.040			1,0



AlphaKite System

www.kometdental.de

Файл / File		SP, rpm	TQ, Ncm	
10.035		500 600	2,0 1,6	
08.040		250-500	2,0	
08.030			1,2	
08.020			0,7	
03.015			0,2	
06.045	1		2,0	
06.040	1		1,5	
06.035	1		1,2	
06.030	1		1,0	
06.025	1		0,7	
06.020	1		0,5	
04.060	2		250-300	3,0
04.055	2		400 500	2,5 2,0
04.050	2	250-500	2,0	
04.045	2		1,8	

Файл / File		SP, rpm	TQ, Ncm
04.040	2	250-500	1,4
04.035	2		0,9
04.030	2		0,7
04.025	2		0,4
04.020	2		0,3
02.035	3		0,7
02.030	3		0,4
02.025	3		0,3
02.020	3		0,2



RaCe
www.fkg.ch

Файл / File		SP, rpm	TQ, Ncm
Pre-RaCe 40/.10		300-600 (500-600)	1,5
Pre-RaCe 35/.08			
Pre-RaCe 40/.06			1,0
Pre-RaCe 30/.06			
40/.06			
35/.06			
30/.06			
25/.06			
20/.06			
15/.06			
10/.06			
50/.04			
40/.04			0,7
35/.04			
30/.04			

Файл / File		SP, rpm	TQ, Ncm
25/.04		300-600 (500-600)	0,5
20/.04			
15/.04			
10/.04			
60/.02			1,0
55/.02			
50/.02			
45/.02			
40/.02			
35/.02			
30/.02			
25/.02			0,5
20/.02			
15/.02			
10/.02			



S-Apex
www.fkg.ch

Файл / File		SP, rpm	TQ, Ncm
ISO 15		500-600	0,5
ISO 20			
ISO 25			
ISO 30			
ISO 35			
ISO 40			
ISO 50			
ISO 60			



iRaCe
www.irace.ch

Файл / File			SP, rpm	TQ, Ncm
R1 (15/.06)			600	1,5
R2 (25/.04)				
R3 (30/.04)				
R1a (20/.02)				
R1b (25/.02)				


















BioRaCe
www.biorace.ch

Файл / File		SP, rpm	TQ, Ncm
BR0 (25/0.08)		500-600	1,0
BR1 (15/0.05)			
BR2 (25/0.04)			
BR3 (25/0.06)			
BR4 (35/0.04)			
BR5 (40/0.04)			
BR4C (35/0.02)			
BR5C (40/0.02)			
BR6 (50/0.04)			
BR7 (60/0.02)			



TWISTED FILES

www.sybronendo.com

Файл / File		SP, rpm	TQ, Ncm	
.12/25		500-600	2,0	
.10/25				
.08/25				
.06/35				
.06/30				
.06/25				
.04/50				
.04/40				
.04/25				



K3™
www.sybronendo.com

Файл / File		SP, rpm	TQ, Ncm
.12/25		350-500	2,8 - 2,0
.10/25			
.08/25			
.06/60			
.06/55			
.06/50			
.06/45			
.06/40			
.06/35			
.06/30			
.06/25			
.06/20			
.06/15			
.04/60			
.04/55			

Файл / File			SP, rpm	TQ, Ncm
.04/50			350-500	3,0 - 2,0
.04/45				
.04/40				
.04/35				
.04/30				
.04/25				
.04/20				
.04/15				
.02/45				
.02/40				
.02/35				
.02/30				
.02/25				
.02/20				
.02/15				



K3™XF
www.sybronendo.com

Файл / File		SP, rpm	TQ, Ncm
.10/25		350-500	2,8 - 2,0
.08/25			
.06/50			
.06/45			
.06/40			
.06/35			
.06/30			
.06/25			
.06/20			
.04/50			
.04/45			
.04/40			
.04/35			
.04/30			
.04/25			
.04/20			



NITI-TEE
www.sendoline.com

Файл / File		SP, rpm	TQ, Ncm
.12/30	Black	300	3,3
.08/30	Blue		
.06/30	Red		3,0
.04/30	Yellow		2,3
.04/25	Yellow		1,2
.04/20	Yellow		0,5
.02/55	Red		
.02/50	Yellow		
.02/45			
.02/40	Black		
.02/35	Green		
.02/30	Blue		
.02/25	Red		
.02/20	Yellow		



S5
www.sendoline.com

Файл / File		SP, rpm	TQ, Ncm
S1 (.08/30)	1	300	3,3
S2 (.06/30)	2		3,0
S3 (.04/30)	3		2,3
S4 (.04/25)	4		1,2
S5 (.04/20)	5		0,5
.04/35	3		2,3
.04/40	3		
.04/45	3		
.04/50	3		
R.12/30	1		
R.08/30	2		3,3
R.06/30	3		3,0
R.04/30	4		2,3
			1,2



HERO 642
www.micro-mega.com

Файл / File		SP, rpm	TQ, Ncm
30/.06	Blue	300-500 600	2,0 1,6
25/.06	Red	300-600	1,5
20/.06	Yellow		1,2
30/.04	Blue		0,8
25/.04	Red		0,5
20/.04	Yellow		0,2
45/.02	White		1,2
40/.02	Black		1,0
35/.02	Green		0,7
30/.02	Blue		0,5
25/.02	Red		0,2
20/.02	Yellow		



HERO SHAPER
www.micro-mega.com

Файл / File		SP, rpm	TQ, Ncm
30/.06	Blue	300-550 600	1,8 1,6
25/.06	Red	300-600	1,3
20/.06	Yellow		1,0
45/.04	White		1,5
40/.04	Black		1,2
35/.04	Green		0,9
30/.04	Blue		0,7
25/.04	Red		0,5
20/.04	Yellow		0,2



ENDOFLARE®
www.micro-mega.com

Файл / File	SP, rpm	TQ, Ncm
ENDOFLARE 25/.12	300-600	2,7 - 1,6






HERO Apical®
www.micro-mega.com

Файл / File	SP, rpm	TQ, Ncm
30/.06	300-600	1,0
30/.08		



R-Endo®






www.micro-mega.com

Файл / File		SP, rpm	TQ, Ncm
Re (25/.12)	_____	300-400	3,0 - 2,5
R1 (25/.08)			2,3
R2 (25/.06)			1,4
R3 (25/.04)			0,7



Revo-S™

www.revo-s.com

Файл / File		SP, rpm	TQ, Ncm
SC1 (25/.06)		250-400	0,8
SC2 (25/.04)			
SU (25/.06)			
AS30 (30/.06)			
AS35 (35/.06)			
AS40 (40/.06)			



CMA System
www.endo-cma.com

Файл / File		SP, rpm	TQ, Ncm
C (25/.08)		300-400	2,5
M (25/.06)			
A1 (20/.04)			
A2 (20/.06)			

DENT
GEO*SOFT*